



Maryland Department of Health Laboratories Administration
INSTRUCTIONS FOR ISOLATE SUBMISSION FOR ARLN AMR REFERENCE TESTING

Isolates must meet one of the following criteria when submitting to the ARLN Regional Laboratory

- Unique gene found in PCR (NDM, OXA-48, IMP, VIM)
- Phenotypic and Genotypic Result Discrepancy
- Outbreak of CRE/CRPA and further characterization warranted

PLEASE DO NOT SEND YOUR CRE/CRPA STATE WIDE SURVEILLANCE SPECIMENS TO THE REGIONAL LABORATORY—EACH STATE JURISDICTIONAL PUBLIC HEALTH LABORATORY SHOULD BE PERFORMING AST AND PCR TESTING FOR THEIR RESPECTIVE SURVEILLANCE SPECIMENS

Notify the MD ARLN Laboratory at mdphl.arln@maryland.gov when sending isolates for Reference Testing – please include the expected receipt date & FedEx tracking number.

FOLLOW THE STEP BY STEP INSTRUCTIONS MENTIONED BELOW:

1. Prepare an agar slant (e.g. TSA, blood agar, etc.) with the isolate to be submitted for AMR reference testing labeled with the appropriate patient identifiers. At least two patient identifiers (such as patient name and date of birth) are required on the specimen under CLIA regulations. Seal the slant or the plate with a piece of parafilm or tape.
2. Go to health.maryland.gov/laboratories/Pages/Home.aspx and click on the Infectious Agents Culture Detection link which will bring up our fillable lab test ordering form which can then be printed after completing or print a scanned copy of the form which can be filled-out manually. Complete the top portion of the form to include all submitting entity and patient demographic information. Include the name of the Healthcare provider who can legally order the test(s) in “Test Request Authorized by” field if these results are to be returned and placed into the patient’s medical record. On the right side, under the section Other Tests for Infectious Agents, or under Reference Microbiology on bottom left side, enter specimen source code in the box and write the name of test as CRE Testing or ARLN Reference testing; also include the name of the isolate species. Attach the provided star shaped large colored “ARLN Reference Test” sticker to the form (See Page 3 for further details). This sticker will help ensure that the specimen/isolate is directed to ARLN Lab for priority reference testing. Without this sticker, the specimen will be processed as a surveillance specimen/isolate, which could be batch tested. More stickers when exhausted can be obtained by contacting the MD ARLN Regional Laboratory at mdphl.arln@maryland.gov. See the attached illustrated directions on page 3 for an example case for completing Form No. 4676. A blank Infectious Agents Culture Detection form No. 4676 is available on page 4. The form must be completed when submitting pre-approved specimens for all AMR reference test requests to the Maryland ARLN Laboratory. Specimens submitted without this form will NOT be accepted for testing. Please ensure that all required core demographic, provider, and patient contact information is completed.
3. Enclose a completed test request form with each specimen/isolate that is submitted for reference testing.
4. Isolates can be transferred within the U.S. as Category B Biological Substances in accordance with Department of Transportation (DoT) Hazardous Materials Regulations (49 CFR Part 171-180). Guidance for packaging samples in accordance with Category B Biological substance requirements can be found in the CDC/NIH Publication Biosafety in Microbiological and BioMedical Laboratories, 5th edition. Additional information about the DoT Hazardous Materials Transport Regulations can be found at <https://www.transportation.gov/pipelines-hazmat>. A supply of pre-printed shipping labels have been provided for your convenience. Appropriately packaged specimens can be shipped directly using the ARLN FEDEX account to the following address:
Maryland Department of Health Laboratories Administration
Attn: ARLN Regional Laboratory
1770 Ashland Ave
Baltimore, Maryland 21205
Email the MD ARLN laboratory at mdphl.arln@maryland.gov about the shipment & the expected receipt date with the FedEx tracking.

STATE LAB
Use Only

Laboratories Administration MDH
1770 Ashland Ave • Baltimore, MD 21205
443-681-3800 <http://health.maryland.gov/laboratories/>
Robert A. Myers, Ph.D., Director



INFECTIOUS AGENTS: CULTURE/DETECTION

HEALTH CARE PROVIDER INFORMATION
HEALTH CARE PROVIDER: TEST FACILITY
Address: 123 N. TEST FACILITY ROAD
City: BALTIMORE County: CITY State: MD Zip Code: 00000
Contact Name: JOE SHMO Phone: 111-111-1111 Fax: 111-111-1112
Test Request Authorized by: DR. TRAB OR SURVEILLANCE/IA

PATIENT INFORMATION
Patient SS # (last 4 digits):
Last name: DOE
First Name: JOHN M.I.
Date of Birth (mm/dd/yyyy): 01 / 01 / 2000
Address: _____ City: _____ State: _____ Zip Code: _____

SEX AND ETHNICITY
Sex: Male Female Transgender M to F Transgender F to M
Ethnicity: Hispanic or Latino Origin? Yes No
Race: American Indian/Alaska Native Asian Black/African American Native Hawaiian/Other Pacific Islander White

MRN/CASE AND DATE INFORMATION
MRN/Case #: _____ DOC #: _____ Outbreak #: _____ Submitter Lab #: _____
Date Collected: 1/1/2017 Time Collected: _____ Onset Date: _____

REASON FOR TEST
Reason for Test: Screening Diagnosis Contact Test of Cure 2-3 Months Post Rx Suspected Carrier Isolate for ID Release
Therapy/Drug Treatment: No Yes Therapy/Drug Type: _____ Therapy/Drug Date: _____

SPECIMEN SOURCE CODE	SPECIMEN SOURCE CODE	SPECIMEN SOURCE CODE
BACTERIOLOGY	MYCOBACTERIOLOGY/AFB/TB	SPECIAL BACTERIOLOGY
Bacterial Culture - Routine	AFB/TB Culture and Smear	Legionella Culture
Add'l Specimen Codes: _____	AFB/TB Referred Isolate for ID	Leptospira
<i>Bordetella pertussis</i>	<i>M. tuberculosis</i> referred Isolate for genotyping	Mycoplasma (Outbreak Investigation Only)
Group A Strep	Nuclear Acid Amplification Test for <i>M. tuberculosis</i> Complex (GeneXpert)	RESTRICTED TESTS Pre-approved submitters only
Group B Strep Screen		
<i>C. difficile</i> Toxin	PARASITOLOGY	<i>Chlamydia trachomatis</i> /GC NAAT
Diphtheria	Blood Parasites: _____	Norovirus** (See comment on reverse)
Foodborne Pathogens (<i>B. cereus</i> , <i>C. perfringens</i> , <i>S. aureus</i>)	Country visited outside US: _____	QuantIFERON
Gonorrhea Culture:	Ova & Parasites	Incubation: Time began: _____ a.m./p.m. Time ended: _____ a.m./p.m.
Incubated? <input type="checkbox"/> Yes <input type="checkbox"/> No	Immigrant? <input type="checkbox"/> Yes <input type="checkbox"/> No	OTHER TESTS FOR INFECTIOUS AGENTS
Hours Incubated: _____	Cryptosporidium	SPEC Test Name: ARLN REFERENCE TESTING
Add'l specimen Codes: _____	Cyclospora/Isospora	Prior arrangements have been made with the following MDH Labs Administration employee: THE MD ARLN TEAM / DR. X
MRSA (rule out)	Microsporidium	SPECIMEN SOURCE CODE PLACE CODE IN BOX NEXT TO TEST
VRE (rule out)	Pinworm	B Blood SP Sputum BW Bronchial Washing T Throat CSF Cerebrospinal Fluid URE Urethra CX Cervix/Endocervix UFV Urine (1 st Void) E Eye UCC Urine (Clean Catch) F Feces V Vagina N Nasopharynx/Nasal W Wound P Penis Other: _____ R Rectum _____ SPECIMEN
ENTERIC INFECTIONS	VIRUS/CHLAMYDIA	
Campylobacter	Adenovirus*	
<i>E. coli</i> 0157 typing/Shiga toxins	<i>Chlamydia trachomatis</i> culture	
Enteric Culture - Routine (Salmonella, Shigella, <i>E. coli</i> 0157, Campylobacter)	Cytomegalovirus	
Salmonella typing	Enterovirus (including coxsackievirus)	
Shigella typing	Influenza	
<i>Vibrio</i>	Respiratory Syncytial Virus	
Yersinia	Parainfluenza (1, 2, 3)	
REFERENCE MICROBIOLOGY	Varicella (VZV)	
ABC's (BIDS) # _____	Organism: _____	
Organism: _____	Comments: _____	
Bacterial Referred Culture for ID		
Specify: ARLN REFERENCE TESTING		

OUR PRINT REQUIRED INFORMATION
PLACE LABELS ON BOTH SIDES

Must complete the submitting entity information (This is where reports will be sent) & include name of practitioner requesting the test.

Complete the patient demographic information such as sex, MRN/Case#, Date collected, time collected, outbreak # and submitter Lab # if available

Patient's first & last names must be on the specimen container and exactly match the test request form

Write "ARLN AMR reference testing" to request testing and enter source code of the specimen/isolate

Check the box and indicate test name by writing CRE Testing or ARLN AMR on the line. Enter Species ID if available.

Using this specimen source code, enter the code in the box next to the test name. If other, please identify the source

Include any comments here

Affix the color coded label



STATE LAB
Use Only

Laboratories Administration MDH
1770 Ashland Ave • Baltimore, MD 21205
443-681-3800 <http://health.maryland.gov/laboratories/>
Robert A. Myers, Ph.D., Director



INFECTIOUS AGENTS: CULTURE/DETECTION

TYPE OR PRINT REQUIRED INFORMATION OR PLACE LABELS ON BOTH COPIES	<input type="checkbox"/> EH <input type="checkbox"/> FP <input type="checkbox"/> MTY/PN <input type="checkbox"/> NOD <input type="checkbox"/> STD <input type="checkbox"/> TB <input type="checkbox"/> CD <input type="checkbox"/> COR		Patient SS # (last 4 digits):		
	Health Care Provider		Last name <input type="checkbox"/> SR <input type="checkbox"/> JR <input type="checkbox"/> Other:		
	Address		First Name M.I.		
	City	County	Date of Birth (mm/dd/yyyy) / /		
	State	Zip Code	Address		
	Contact Name:		City County		
	Phone #	Fax #	State Zip Code		
	Test Request Authorized by:				
	Sex: <input type="checkbox"/> Male <input type="checkbox"/> Female <input type="checkbox"/> Transgender M to F <input type="checkbox"/> Transgender F to M		Ethnicity: Hispanic or Latino Origin? <input type="checkbox"/> Yes <input type="checkbox"/> No		
	Race: <input type="checkbox"/> American Indian/Alaska Native <input type="checkbox"/> Asian <input type="checkbox"/> Black/African American <input type="checkbox"/> Native Hawaiian/Other Pacific Islander <input type="checkbox"/> White				
MRN/Case #	DOC #	Outbreak #	Submitter Lab #		
Date Collected:	Time Collected:	<input type="checkbox"/> a.m. <input type="checkbox"/> p.m.	Onset Date: / /		
Reason for Test: <input type="checkbox"/> Screening <input type="checkbox"/> Diagnosis <input type="checkbox"/> Contact <input type="checkbox"/> Test of Cure <input type="checkbox"/> 2-3 Months Post Rx <input type="checkbox"/> Suspected Carrier <input type="checkbox"/> Isolate for ID <input type="checkbox"/> Release					
Therapy/Drug Treatment: <input type="checkbox"/> No <input type="checkbox"/> Yes		Therapy/Drug Date: / /			
↓ SPECIMEN SOURCE CODE		↓ SPECIMEN SOURCE CODE		↓ SPECIMEN SOURCE CODE	
↓ BACTERIOLOGY		↓ MYCOBACTERIOLOGY/AFB/TB		↓ SPECIAL BACTERIOLOGY	
Bacterial Culture - Routine		AFB/TB Culture and Smear		Legionella Culture	
Add'l Specimen Codes: _____		AFB/TB Referred Isolate for ID		Leptospira	
<i>Bordetella pertussis</i>		<i>M. tuberculosis</i> referred Isolate for genotyping		Mycoplasma (Outbreak Investigation Only)	
Group A Strep		Nuclear Acid Amplification Test for <i>M. tuberculosis</i> Complex (GeneXpert)		RESTRICTED TESTS Pre-approved submitters only	
Group B Strep Screen					
<i>C. difficile</i> Toxin		PARASITOLOGY		<i>Chlamydia trachomatis</i> /GC NAAT	
Diphtheria		Blood Parasites:		Norovirus** (See comment on reverse)	
Foodborne Pathogens		Country visited outside US:		QuantiFERON	
<i>(B. cereus, C. perfringens, S. aureus)</i>		Ova & Parasites		Incubation: Time began: ____ a.m./p.m.	
Gonorrhea Culture:		Immigrant? <input type="checkbox"/> Yes <input type="checkbox"/> No		Time ended: ____ a.m./p.m.	
Incubated? Yes ___ No		Cryptosporidium		OTHER TESTS FOR INFECTIOUS AGENTS	
Hours Incubated: _____		Cyclospora/Isospora			
Add'l specimen Codes: _____		Microsporidium		Test Name: _____	
MRSA (rule out)		Pinworm		Prior arrangements have been made with the following MDH Labs Administration employee: _____	
VRE (rule out)		VIRUS/CHLAMYDIA			
ENTERIC INFECTIONS		Adenovirus*			
Campylobacter		<i>Chlamydia trachomatis</i> culture			
<i>E. coli</i> 0157 typing/Shiga toxins		Cytomegalovirus (CMV)			
Enteric Culture - Routine		Enterovirus (Includes Echo & Coxsackie)		SPECIMEN SOURCE CODE PLACE CODE IN BOX NEXT TO TEST	
<i>(Salmonella, Shigella, E. coli 0157, Campylobacter)</i>		Herpes Simplex Virus (Types 1 & 2)			
Salmonella typing		Influenza (Types A & B)* Rapid Flu Test:			
Shigella typing		Type: _____			
<i>Vibrio</i>		Result: <input type="checkbox"/> Negative <input type="checkbox"/> Positive			
Yersinia		Patient admitted to hospital? <input type="checkbox"/> Yes <input type="checkbox"/> No			
REFERENCE MICROBIOLOGY		Parainfluenza (Types 1, 2 & 3)*			
ABC's (BIDS) # _____		Varicella (VZV)			
Organism: _____		*MAY INCLUDE RESPIRATORY SCREENING PANEL			
Bacteria Referred Culture for ID		Comments:			
Specify: _____					
				B Blood SP Sputum	
				BW Bronchial Washing T Throat	
				CSF Cerebrospinal Fluid URE Urethra	
				CX Cervix/Endocervix UFV Urine (1 st Void)	
				E Eye UCC Urine (Clean Catch)	
				F Feces V Vagina	
				N Nasopharynx/Nasal W Wound	
				P Penis O Other: _____	
				R Rectum _____	

DHMH 4676 Revised 05/17

Client