

Isolation and Identification of Legionella at the MD DHMH Laboratories Administration

Thomas Lawson, M.S. and Dr. Maria Paz Carlos
Division of Virology and Immunology
**State of Maryland Department of Health and Mental
Hygiene Laboratories Administration**



Clinical Tests for Legionella

Performed at the MD DHMH Laboratories

- **Culture/ Direct Fluorescent Antibody (i.e., sputum, lung tissue)**
- **Serology IFA**
- **Urine Antigen EIA/ Rapid**

Additional informational found at MD DHMH Guide to Public Health Lab Services: <http://dhmh.maryland.gov/laboratories/docs/guide.pdf>

Clinical Test Request Submission Form for Legionella



Laboratories Administration MD DHMH
 201 W. Preston St. • Baltimore, MD 21201
 P.O. Box 2355 • Baltimore, MD 21203-2355
 410-767-6100 www.dhmm.state.md.us/labs
 Robert A. Myers, Ph.D., Director

STATE LAB
Use Only

INFECTIOUS AGENTS: CULTURE/DETECTION

TYPE OR PRINT REQUIRED INFORMATION OR PLACE LABELS ON ALL FOUR COPIES	DEH <input type="checkbox"/> DFP <input type="checkbox"/> DMTY/PN <input type="checkbox"/> DNOD <input type="checkbox"/> DSTD <input type="checkbox"/> DTB <input type="checkbox"/> OCD <input type="checkbox"/> DCOR <input type="checkbox"/>		Patient SS# (last 4 digits):	
	Submitter		Last Name <input type="checkbox"/> SR <input type="checkbox"/> JR <input type="checkbox"/> Other <input type="checkbox"/>	
	Address		First Name M.I. Maiden:	
	City County		Date of Birth (mm/dd/yyyy) / /	
	State Zip Code		Address	
	Contact Name		City County	
	Phone # Fax #		State Zip Code	
	Ethnicity: <input type="checkbox"/> Hispanic/Latino <input type="checkbox"/> Not Hispanic/Latino <input type="checkbox"/> Unknown		Sex: <input type="checkbox"/> Male <input type="checkbox"/> Female <input type="checkbox"/> Transgender	
	Race: <input type="checkbox"/> White <input type="checkbox"/> African American <input type="checkbox"/> Asian/Pacific <input type="checkbox"/> American Indian/Alaska Native <input type="checkbox"/> Multiracial <input type="checkbox"/> Not Specified <input type="checkbox"/> Other		Case # DOC# Outbreak # Submitter Lab#	
	Collect Date: Collect Time: <input type="checkbox"/> am <input type="checkbox"/> pm Onset Date:		Name of Test Date <input type="checkbox"/> 1 st <input type="checkbox"/> 2 nd <input type="checkbox"/> 3 rd State Lab Number:	
Reason for Test: <input type="checkbox"/> Screening <input type="checkbox"/> Diagnosis <input type="checkbox"/> Contact <input type="checkbox"/> Test of Cure <input type="checkbox"/> 2-3 Months Post Rx <input type="checkbox"/> Suspected Carrier <input type="checkbox"/> Isolate for ID <input type="checkbox"/> Release		Therapy/Drug Treatment: <input type="checkbox"/> No <input type="checkbox"/> Yes Therapy/Drug Type: _____		

SPECIMEN CODE	SPECIMEN CODE	SPECIMEN CODE
BACTERIOLOGY	SPECIAL BACTERIOLOGY	RESTRICTED TESTS
Bacterial Culture - Routine	Legionella Culture	Pre-approved submitters only
Additional specimen codes:	Leptospira	Chlamydia trachomatis/GC NAAT
<i>Bordetella pertussis</i>	Mycoplasma	Chlamydia trachomatis only/NAAT
Group A Strep	MYCOBACTERIOLOGY/AFB/TB	Norovirus ** (see comment on back)
Group B Strep Screen	AFB/TB Culture and Smear	OTHER TESTS FOR INFECTIOUS AGENTS
<i>C. difficile</i> Toxin	AFB/TB Referred Culture for ID	Test name: _____
Diphtheria	<i>M. tuberculosis</i> Referred Culture for Genotyping	Prior arrangements have been made with the following DHMH Laboratories Administration employee: _____
Foodborne Pathogens (<i>B. cereus</i> , <i>C. perfringens</i> , <i>S. aureus</i>)	Nucleic Acid Amplification Test for <i>M. tuberculosis</i> Complex (MTD)	
Gonorrhea Culture: Incubated? <input type="checkbox"/> Yes <input type="checkbox"/> No	PARASITOLOGY	
Hrs. incubated: _____ Add'l specimen codes: _____	Blood Parasites:	
MRSA (rule out)	Country visited outside US: _____	SPECIMEN CODE: PLACE CODE IN BOX NEXT TO TEST
VRE (rule out)	Ova & Parasites: Immigrant? <input type="checkbox"/> Yes <input type="checkbox"/> No	B Blood
ENTERIC INFECTIONS	Cryptosporidium	BW Bronchial Washing
Campylobacter	Cyclospora/Isospora	CSF Cerebrospinal Fluid
<i>E. coli</i> O157 typing	Microsporidium	CX Cervix/Endocervix
Enteric Culture - Routine (Salmonella, Shigella, <i>E. coli</i> O157, Campylobacter)	Pinworm	E Eye
Salmonella typing	VIRUS/CHLAMYDIA	F Feces
Shigella typing	Adenovirus*	N Nasopharynx/Nasal
<i>V. parahaemolyticus</i>	Arbovirus Panel (WNV, EEEV, SLEV)	P Penis
Yersinia	Chlamydia trachomatis	R Rectum
REFERENCE MICROBIOLOGY	Cytomegalovirus (CMV)	SP Sputum
ABC'S (BIDS) #	Enterovirus (Inc. Echo & Coxsackie)	T Throat
Organism: _____	Herpes Simplex Virus (Types 1 & 2)	URE Urethra
Bacteria Referred Culture for ID	Influenza (Types A & B)*	UFV Urine (First Void)
Specify: _____	Parainfluenza (Types 1, 2 & 3)*	VCC Urine (Clean Catch)
	Respiratory Syncytial Virus (RSV)*	W Vagina
	Varicella (VZV)	O Wound
		O Other:
	*MAY INCLUDE RESPIRATORY SCREENING PANEL.	
	Comments: _____	



DHMH 4676 Revised 09/11

ORIGINAL



Laboratories Administration MD DHMH
 201 W. Preston St. • Baltimore, MD 21201
 P.O. Box 2355 • Baltimore, MD 21203-2355
 410-767-6100 www.dhmm.state.md.us/labs
 Robert A. Myers, Ph.D., Director

STATE LAB
Use Only

SEROLOGICAL TESTING

TYPE OR PRINT REQUIRED INFORMATION OR PLACE LABELS ON ALL FOUR COPIES	DEH <input type="checkbox"/> DFP <input type="checkbox"/> DMTY/PN <input type="checkbox"/> DNOD <input type="checkbox"/> DSTD <input type="checkbox"/> DTB <input type="checkbox"/> OCD <input type="checkbox"/> DCOR <input type="checkbox"/>	Patient SS# (last 4 digits):	
	Submitter	Last Name <input type="checkbox"/> SR <input type="checkbox"/> JR <input type="checkbox"/> Other <input type="checkbox"/>	
	Address	First Name M.I. Maiden:	
	City County	Date of Birth (mm/dd/yyyy) / /	
	State Zip Code	Address	
	Contact Name	City County	
	Phone # Fax #	State Zip Code	
	Ethnicity: <input type="checkbox"/> Hispanic/Latino <input type="checkbox"/> Not Hispanic/Latino <input type="checkbox"/> Unknown	Sex: <input type="checkbox"/> Male <input type="checkbox"/> Female <input type="checkbox"/> Transgender	
	Race: <input type="checkbox"/> White <input type="checkbox"/> African American <input type="checkbox"/> Asian/Pacific <input type="checkbox"/> American Indian/Alaska Native <input type="checkbox"/> Multiracial <input type="checkbox"/> Not Specified <input type="checkbox"/> Other	Case # DOC# Outbreak # Submitter Lab#	
	Collect Date: Collect Time: <input type="checkbox"/> am <input type="checkbox"/> pm *Rabies Vaccination Dates	Name of Test Date <input type="checkbox"/> 1 st <input type="checkbox"/> 2 nd <input type="checkbox"/> 3 rd State Lab Number:	
Previous Test Done? <input type="checkbox"/> No <input type="checkbox"/> Yes	Name of Test Date <input type="checkbox"/> 1 st <input type="checkbox"/> 2 nd <input type="checkbox"/> 3 rd State Lab Number:		
Onset Date: Exposure Date: Clinical Illness:			

SEROLOGY: SERUM (1ml/test) or WHOLE BLOOD (5ml) REQUIRED	LAVENDER TOP TUBE REQUIRED
Arbovirus / West Nile Virus Panel	Herpes Simplex Virus (HSV) Types 1&2
(Serum or CSF)	Legionella
Provide dates of onset & collection (see above)	Leptospira <input type="checkbox"/> Yes <input type="checkbox"/> No
Required information, check all that apply:	Lyme Disease <input type="checkbox"/> Yes <input type="checkbox"/> No
DIAGNOSIS: <input type="checkbox"/> aseptic meningitis <input type="checkbox"/> encephalitis	MMEV Immunity Screen: <input type="checkbox"/> Measles (Rubeola), Mumps, Rubella, Varicella (Chickenpox) IgG Ab only
<input type="checkbox"/> fever <input type="checkbox"/> other _____	Mononucleosis - Infectious
SYMPTOMS: <input type="checkbox"/> headache <input type="checkbox"/> fever <input type="checkbox"/> stiff neck	Mumps Immunity Screen
<input type="checkbox"/> altered mental status <input type="checkbox"/> muscle weakness	Mycoplasma
<input type="checkbox"/> rash <input type="checkbox"/> other _____	Rocky Mountain Spotted Fever (RMSF)
ILLNESS FATAL? <input type="checkbox"/> Yes <input type="checkbox"/> No	Rabies (RFFIT) (*List vaccination dates above)
TRAVEL HISTORY (dates and places)	Rubella Immunity Screen
IMMUNIZATIONS: Yellow fever? <input type="checkbox"/> Yes <input type="checkbox"/> No	Rubella (Measles) Immunity Screen
Flavivirus? <input type="checkbox"/> Yes <input type="checkbox"/> No	Syphilis - Previously treated? <input type="checkbox"/> Yes <input type="checkbox"/> No
IMMUNOCOMPROMISED? <input type="checkbox"/> Yes <input type="checkbox"/> No	Toxoplasma
Aspergillus	Tularemia
Brucella	Varicella Immunity Screen
Chlamydia (group antigen IgG)	VDRL (CSF only)
Cryptococcal antigen)	CDC Test request
Cytomegalovirus (CMV)	Call lab (410-767-6162) prior to submitting
<i>E. histolytica</i>	Test requested: _____
Ehrlichia	Other test request: _____
Epstein-Barr Virus (EBV)	
Hepatitis A Screen (IgM Ab only, acute infection)	
Call lab (410-767-6169) prior to submitting	
Hepatitis B Screen (HBs antigen only)	
Prenatal patient? <input type="checkbox"/> Yes <input type="checkbox"/> No	
Hepatitis B Panel: (HBsAg, HBsAb)	
Reflex Testing: HBsAb Neg HB Total Core	
HBsAg Pos: HB Core IgM, HBsAg, HBeAb	
Hepatitis B post vaccine	
Hepatitis C Screen (HCV Ab only)	
	Prior arrangements have been made with the following DHMH Labs Administration employee: _____



ORIGINAL





DHMH 4677 Revised 11/11

Environmental Tests for Legionella Performed at the MD DHMH Laboratories

Isolation by Culture

Identification by Direct Fluorescent Antibody (DFA) Test

MD DHMH Laboratories Administration – **Certified by CDC
ELITE (Environmental Legionella Isolation Techniques
Evaluation)**



Legionella Possible Case Scenario in Maryland: The Call

- Lab receive call from State Epi
- Preapproval
- Pre-remediation
- Post-remediation
- Coordination with Epi, HDs
- Collection Sites in Facility
- Collection Kit
- Specimen Submission Policy



Sample Collection Kit

- Sterile 1 L bottles
- Sodium thiosulfate
- Swabs
- Sterile 5ml tubes
- Sterile water
- DHMH 4676
Culture Detection
Submission Forms

➤ Storage Condition for collection kit: 4°C



Bulk Water Sample Collection

- **Bulk water: Minimum of 1 liter**
- **Add 0.5 ml of 0.1N sodium thiosulfate to each 1L water bottle.**




Swab Sample Collection

- Swabs of faucet aerators and shower heads should be taken before water samples from these sites.
- The sample should be taken with the aerator or shower head removed if possible.
- Submerge in 3-5 ml of water taken at the same time to prevent drying during transport to the laboratory.



Specimen Packaging

- **Collected Samples**
Ensure that the water bottles and tubes are tightly closed.
- **Completed Submission Forms (DHMH 4676 Culture Detection)**



Sample Handling and Delivery


Sample delivery MUST be coordinated with the laboratories prior to submission. Please call 410-767-6162 (Supervisor: Thomas Lawson).

If samples cannot be delivered directly to the lab after collection, store at 4°C.

Samples must be received within 48 hours after collection.

MD DHMH Laboratories Administration

<http://dhmh.state.md.us/laboratories>



Samples Received at the Lab

Samples Acceptability Criteria



- Collected samples
- Completed DMMH 4676 Culture Detection Submission Forms
- Printed/e-mail Sample Log

Laboratory Information Management System

STARLIMS

STARLIMS MASTER V. 2.006 (MSSQL) , User Name - LAWSONT, Language - ENG

File Design Help

Console

- Patient Demographics
- Animal Data
- Food Data
- BT Data
- Environmental Data
- Micro Reports
- Bench Log Reports
- Final Log Report
- Micro Final Log Report
- Print Untested Samples
- Print Rejection Report
- Bio T Manual Labels Print
- Print Unreleased Samples
- Print Released Line List
- Rabies Manual Labels Print
- Re-Print Barcode Labels
- TB Manual Labels Print
- Registration/Reporting
- Critical Link Statistics
- Syphilis Reports
- Childhood Screening
- Norovirus Reports
- Molecular Reports
- Chlamydia
- Microbiology
- EDCP Reports
- Virology
- Environmental
- Edit**
- Clinical Samples
- View Samples
- Environmental Samples
- Released Clinical Sample
- Released Env. Samples
- Edit run
- View Audit**
- Audit HL7
- Audit Records
- DII Queue

Select

External Patient ID	<input type="text"/>	Internal Sample ID	<input type="text"/>	Accession#	<input type="text"/>
UI Number	<input type="text"/>	CTR Number	<input type="text"/>	DOC ID Number	<input type="text"/>
SS Number	<input type="text"/>	Birthdate:	<input type="checkbox"/>	Tuesday , May 08, 2012	<input type="text"/>
First Name	<input type="text"/>	Last Name	<input type="text"/>		
Animal Name	<input type="text"/>		Run#:	<input type="text"/>	
CQ#:	<input type="text"/>		Sub Lab #:	<input type="text"/>	
Owner First Name	<input type="text"/>		Owner Last Name	<input type="text"/>	
EID	<input type="text"/>	Rabies #	<input type="text"/>	Bite#	<input type="text"/>
Outbreak #	<input type="text"/>				
Conditions/Tests	<input type="text"/>				
Submitters	<input type="text"/>				
Date Logged	<input type="checkbox"/>	From	Tuesday , May 08, 2012	To	Tuesday , May 08, 2012
Date Received	<input type="checkbox"/>	From	Tuesday , May 08, 2012	To	Tuesday , May 08, 2012

Ok Cancel

Sample Processing


Bulk Water

Filtration

- 1L bulk water
- 0.2 μ M polycarbonate filter

Swab

- Direct plating



Inoculation



Inoculate into 3 types of culture plates.

- 1. PAV (Polymyxin -B Anisomycin Vancomycin)**
- 2. EBCYE (Environmental Buffered Charcoal Yeast Extract)**
- 3. BCYE (Buffered Charcoal Yeast Extract)**

Incubation



37°C Incubation

2.5% CO₂

Plates checked daily after Day 2


Examination of Culture Plates

Collection of Cultures for Legionellae


**Colonies convex
round with entire
edges**

**Center of the colony
usually bright white
with textured edges**


**Once identified,
individual colonies
are collected with
loop**



Colony Screen: DFA (Direct Fluorescent Antibody)



Cell Suspensions




Staining



Antibodies and Antigen Controls




Initial Screen



***Legionella pneumophila*
serogroups 1-14**

***Legionella non-pneumophila*
spp. (b-m)**

DFA Fluorescent Microscopy




Reactive

Legionella pneumophila Serogroups 1-14








Continue with *L. pneumophila* subculture,
serotyping, blood agar plate (BAP)




Non-Reactive

Legionella non-pneumophila species (b-m)

Legionella pneumophila Serotyping




Continue with blood agar plate inoculation



Non-Reactive

Blood Agar Plate (BAP) Confirmation



Growth on BAP

**Absence of
Legionella**

24 hour
incubation

No Growth on BAP

**Presence of
Legionella**

**Calculate CFU/ml on
original plates**



State of Maryland
Department of Health and Mental Hygiene
LABORATORIES ADMINISTRATION

Central Laboratory
201 West Preston Street
Baltimore, MD 21201
Robert Myers, Ph.D., Director
www.dhmd.state.md.us/labs

SAMPLE REPORT

201 W PRESTON STREET
3RD FLOOR, O'CONNOR BUILDING
BALTIMORE, MD 21201

Specimen Number:
Patient Name:
Patient Id:
Birth Date:
Sex:
SSN:
Patient Address:

Internal Lab Number:
Outbreak Number:
Specimen Source:
Date Collected:
Date/Time Received:
Submitter Lab #

Comments: LEGIONELLA OUTBREAK INVESTIGATION - SWAB

<u>Test Name</u>	<u>Results</u>	<u>Date Reported</u>
Legionella Culture Isolation	L. pneumophila serogroup 4 isolated after 10 days incubation.	
Legionella Culture Isolation	L. pneumophila serogroup 5 isolated after 10 days incubation.	
Legionella Culture Isolation	L. pneumophila serogroup 6 isolated after 10 days incubation.	
Legionella Culture Isolation	No L. pneumophila serogroup 1 to 6 or L. micdadei isolated after 10 days incubation.	
Legionella Culture Isolation	L. pneumophila isolated after 10 days incubation.	
Legionella Culture Isolation	Legionella species isolated after 10 days incubation.	

Comments:

CFU/ML = 3.39

Conclusions


- **Proper specimen collection and submission are essential for quality assurance testing.**
- **Clinical and Environmental Case Linkage**
- **Strong coordination and collaboration with Epi and Local Health Departments are critical.**

Future Directions

Collaboration with Molecular on developing a PCR assay for Lp1.

Legionella Laboratory Website

MyLIMS Legionella Roll-Out



Acknowledgements

Maryland DHMH

David Blythe

Maria Carlos

Gonzalo Crujeiras

Alvina Chu

Ni Lui

Robert Myers

Rene Najera

Vanlila Patel

Maria Said

Ruth Thompson

Leena Trivedi

Lucy Wilson

